

PLASMID DNA EXTRACTION MINI KIT

Description

The DNA Plasmid Extraction Kit provides a rapid, phenol-free method for the extraction of high-purity plasmid DNA from bacterial cultures such as *E. coli*, which bacteria is pellet, lysed, and then neutralized. The extracted DNA can be used in a variety of applications such as PCR, cloning, sequencing, in vitro transcription, and labeling. Also, as a column-type tube is utilized in the purification process, extraction is carried out in three simple steps of binding / washing / elution. Once bound, the DNA is washed and then eluted from the column, ready for use.

Features

- For high yields of plasmid DNA-up to 30µg from 1~5ml overnight cultures.
- Effective purification of DNA fragments ranging from 100bp to 12+kb.
- No need for messy resin slurries, extracting with phenol, or concentrating via alcohol precipitation.
- Superior purity-DNA yields quality sequence data using automated or manual methods.
- Optimized buffers are included for maximum DNA purity and yield.
- Versatile protocol-works with all neutral gel buffers and both conventional and low-melting agarose gel.

Format

Spin Columns

Operation

Centrifuge/vacuum

Binding Capacity

Up to 30 µg

Expectant Yield

20-30 µg for high-copy plasmid
3-10 µg for low-copy plasmid

Time Required


Under 30 minutes

Applications

- Fluorescent or radioactive Sequencing
- Ligation
- Restriction enzyme digestion
- Ligation and Transformation
- Library screening

Time Required


Under 30 minutes



For Research Use Only

DNA Sequencing Analysis

Sequencing of plasmid purified with Plasmid DNA Extraction Mini Kit. The sequence was analyzed with BigDye Terminator chemistry on ABI 3700




Ordering Information

Cat. No.	Product Name	Size	Kit Components	Store at
DE-005	Plasmid DNA Extraction Kit	100 preps.	<ol style="list-style-type: none"> 1. FSPD 1 Buffer 2. FSPD 2 Buffer 3. FSPD3 Buffer 4. W1 Buffer 5. Wash Buffer (conc.) 6. Elution Buffer 7. RNase A (50 mg/ml) 8. FSPD Columns 9. 2 ml Collection Tube 	FSPD 1 Buffer with RNaseA at 4 °C, other components should at room temperature (15 -25 °C).
DE-006	Plasmid DNA Extraction Kit	300 preps.	<ol style="list-style-type: none"> 1. FSPD 1 Buffer 2. FSPD 2 Buffer 3. FSPD3 Buffer 4. W1 Buffer 5. Wash Buffer (conc.) 6. Elution Buffer 7. RNase A (50 mg/ml) 8. FSPD Columns 9. 2 ml Collection Tube 	FSPD 1 Buffer with RNaseA at 4 °C, other components should at room temperature (15 -25 °C).

For Research Use Only

M1 1 2 3 4 5 6 M2



Agarose Gel Analysis Of Plasmids Purified

Agarose gel analysis of 5 plasmid and cosmid extracted by the Plasmid DNA Extraction Mini Kit.

M: 1kb DNA ladder (100 bp-12,000 bp)

1: pBluescript II

2: pUC 18

3: pBR 322

4: pGem with 1 Kb insert

5: pET 43.1 with 2.7 Kb insert

6: cosmid (50 Kb)

M1: 1 Kb DNA Ladder

M2: Lambda-HindIII

Procedure

The modified alkaline lysis method and RNase treatment are used for creating cleared cell lysate with minimal genomic DNA and RNA contaminants. In the presence of a chaotropic salt, the plasmid DNA in the lysate binds to the glass fiber matrix in the spin column. The purified plasmid DNA is eluted by a low salt elution buffer or water. The procedure doesn't require DNA phenol extraction or alcohol precipitation.

Storage Conditions

Plasmid DNA Extraction Kit can be stored at room temperature (15-25°C). After adding RNase A, FSPD Buffer1, should be stored at 4°C and is stable for six months. Other buffers and columns can be stored dry for up to 1 year at room temperature (15-25°C).

Distributed by:



**SOCIETÀ
ITALIANA
CHIMICI**
DIVISIONE SCIENTIFICA r.l.

Via Rio nell'Elba 140 - 00138 Rome - ITALY
Phone: +39 06 8818936 / 6800211 Fax: +39 06 8815319
Web-Site: www.sichim.com E-mail: info@sichim.com

For Research Use Only

FISHER MOLECULAR BIOLOGY

36 Terry Drive
Trevose, PA 19048 - USA

Introduction

Plasmid Extraction Mini Kit is an excellent tool offering a speed and economic method to purify plasmid DNA from bacteria cultures. This technology is based on binding DNA to silica-based membranes in chaotropic salts, washing DNA with ethanol-contained Wash Buffer. Compare with other harmful and time-consuming procedure, such as phenol/ chloroform extraction and ethanol precipitation, Plasmid extraction kit shortens the handling time to about 25 minutes. The high quality plasmid DNA can be used directly for the downstream application.

Specification

Sampling: 1~5 ml overnight culture

Plasmid Size: < 12Kb

Yield: 20~30 µg of high-copy plasmid

Handling time: about 25 min

Kit Contents:	DE-005	DE-006
FSPD1 Buffer	30 ml	90 ml
FSPD2 Buffer	30 ml	90 ml
FSPD3 Buffer	40 ml	120 ml
Wash Buffer (concentration)**	20 ml	50 ml
W1 Buffer (concentration)*	35 ml	98 ml
Elution Buffer	15 ml	35 ml
RNase A (50mg/ml)	60 µl	180 µl
FSPD Column	100 pcs	300 pcs
Collection Tube	100 pcs	300 pcs
User Manual	1	1

* Add 13 ml/ 36 ml ethanol (96 ~ 100%) to W1 Buffer when first open.


** Add 80 ml/ 200 ml ethanol (96 ~ 100%) to Wash Buffer when first open.

For Research Use Only

Important Notes:

1. Buffer provided in this kit contain irritants. Wear gloves and lab coat when handling these buffer.
2. Brief spin RNase A tube to remove drops from the inside of the lid. Add 1 ml of FSPD1 Buffer into RNase A tube and mix well. Transfer the mixture into FSPD1 Buffer bottle and store at 4 ° C.
3. Check FSPD2 Buffer before use. Warm FSPD2 Buffer at 55 ° C for 10 minutes if any precipitation formed. Don't shake FSPD2 Buffer vigorously.
4. To avoid acidification of FSPD2 Buffer from CO₂ in the air, close the bottle immediately after use.
5. For DE-005, add 13 ml ethanol (96~100%) to W1 Buffer when first open. For DE-006, add 36 ml ethanol (96~100%) to W1 Buffer when first open.
6. For DE-005, add 80 ml ethanol (96~100%) to Wash Buffer when first open. For DE-006, add 200 ml ethanol (96~100%) to Wash Buffer when first open.
7. All centrifuge steps are done at full speed (14,000 rpm or 10,000 x g) in a microcentrifuge.

Brief Procedure



For Research Use Only

• **General Protocol:**

1. Transfer 1-5 ml of well-grown bacteria culture to a microcentrifuge tube (not provided).
2. Descend the bacteria by centrifuging for 1-2 min and discard the supernatant completely.
3. Add 250 µl of FSPD1 Buffer to the pellet and resuspend the cells completely by pipetting.
 - **Make sure that RNase A has been added into FSPD1 Buffer when first open.**
 - **No cell pellet should be visible after resuspension of the cells.**
4. Add 250 µl of FSPD2 Buffer and gently invert the tube 5 times to lyse the cells and incubate at room temperature for 2 min.
 - **Do not vortex, vortex may shear genomic DNA. If necessary, continue inverting the tube until the lysate become clear.**
 - **Do not proceed this step over 5 min.**
5. Add 350 µl of FSPD3 Buffer and invert the tube 5 times immediately but gently.
 - **Invert immediately after add FSPD3 Buffer will avoid asymmetric precipitation.**
6. Centrifuge for 10 min. During centrifuging, place a FSPD Column in a Collection Tube.
7. Transfer the supernatant carefully to FSPD Column. Centrifuge for 1 min then discard the flow-through.
 - **Do not transfer any white pellet into the column.**
8. Add 400 µl of W1 Buffer to FSPD Column. Centrifuge for 1 min then discard the flow-through.
 - **Make sure that ethanol (96-100 %) has been added into W1 Buffer when first open.**
9. Add 750 µl of Wash Buffer to FSPD Column. Centrifuge for 1 min then discard the flow-through.
 - **Make sure that ethanol (96-100 %) has been added into Wash Buffer when first open.**
10. Centrifuge for an additional 5 min to dry the column.
 - **Important step ! This step will remove the residual liquid completely that will inhibit subsequent enzymatic reaction.**
11. Place FSPD Column to a new 1.5 ml microcentrifuge tube (not provided).
12. Add 50 µl ~ 100 µl of Elution Buffer or ddH₂O to the membrane center of FSPD Column. Stand the column for 1 min.
 - **Important step ! For effective elution, make sure that the elution solution is dispensed on the membrane center and is absorbed completely.**

- **Important : Do not elute the DNA using less than suggested volume (50µl). It will lower the final yield.**
13. Centrifuge for 1 min to elute plasmid DNA.
14. Store plasmid DNA at 4 °C or -20 °C.

Distributed by:



Via Rio nell'Elba 140 - 00138 Rome - ITALY
Phone: +39 06 8818936 / 8800211 Fax: +39 06 8815319
Web-Site: www.sichim.com E-mail: info@sichim.com

For Research Use Only

FISHER MOLECULAR BIOLOGY

36 Terry Drive
Trevose, PA 19048 - USA

• Troubleshooting

Low yield

Bacterial cells were not lysed completely

- **Too many bacterial cells were used (OD600 > 10). Separate the bacterial culture into multiple tubes.**
- **After FSPD3 Buffer addition, break up the precipitate by inverting to ensure higher yield.**

Overgrown of bacterial cells

- **Incubation time should not longer than 16 hours.**

Bacterial cells were insufficient

- **Ensure that bacterial cells have grown to an expected amount (OD600 > 1) after incubation under suitable shaking modes.**

Incorrect DNA Elution Step

- **Ensure that Elution Buffer was added and absorbed to the center of FSPD Column Martix.**

Incomplete DNA Elution

- **If size of DNA fragments is larger than 10 kb, use preheated Elution Buffer (60~70°C) on Elution Step to improve the elution efficiency.**

Incorrect Wash Buffer

- **Ensure that Ethanol was added to Wash Buffer prior to use.**
-

Eluted DNA does not perform well

Residual ethanol contamination

- **After Wash Step, dry FSPD Column with additional centrifugation at top speed for 5 minutes or incubation at 60°C for 5 minutes.**
-

Genomic DNA Contaminates

Lysate prepared improperly.

- **Gently invert the tube after adding FSPD2 Buffer. And the incubation time should not longer than 5 minutes.**
 - **Do Not use overgrown bacterial culture.**
-

• Troubleshooting

RNA Contaminates Plasmid DNA

Insufficiency of RNase A activity in FSPD1 Buffer because of long-term storage

- **Prior to using FSPD1 Buffer, ensure that RNase A was added. If RNase A added FSPD1 Buffer is out of date, add additional RNase A into FSPD1 Buffer to a concentration of 50µg/ml then store 4°C.**
 - **Too many bacterial cells were used, reduce sample volume.**
-

Smearing or degrading of Plasmid DNA

Nuclease contamination

- **If used host cells have high nuclease activity (e.g., enA+ strains), perform this Optional Wash Step to remove residuary nuclease.**
 - **After DNA Binding Step, add 400µl of W1 Buffer into FSPD column and column and incubate for 2 minutes at room temperature.**
 - **Centrifuge at full speed (14,000 rpm or 10,000 xg) for 30 seconds.**
 - **Followed using standard Wash Step.**
-

Plasmid DNA is not adequate for enzymatic digestions

Eluted plasmid DNA contains residual ethanol

- **Make sure you have discarded the flow-through after washing with Wash Buffer (Step 9) and centrifuged for an additional 3 minutes (Step 10).**
-

Denatured Plasmid DNA migrate faster than supercoiled form during electrophoresis

Incubation in FSPD2 Buffer is too long

- **Do not incubate longer than 5 minute in FSPD2 Buffer**
-

For Research Use Only