

FS-T-70222 TaqMan Universal Multiplex qPCR Master Mix

Description

TaqMan Universal Multiplex PCR Master Mix is a ready to use master solution for multiple 5' nuclease DNA applications. Supplied at 2X concentration. The mix contains Hot Start DNA Polymerase, Dnase inhibitor, MgCl₂, a passive internal reference based on proprietary ROX dye. It contains buffer enhancements to guarantee performance and reliability in all including multiplex fluorescence quantitative experiments.

- Rox Dye are supplied as a separate vials.

This master mix can be used for DNA detection: Genomic DNA, cDNA, plasmid DNA, and Viral Sequences.

This Master Mix can be used for gene typing and gene multiplex quantitative analysis.

| Description | FS-T-70222 |
|---|-------------|
| 2X TaqMan Universal Multiplex qPCR Universal Master Mix | 2 x 1.25 mL |
| 50x Rox Dye I | 100 ul |
| 50x Rox Dye II | 100 ul |

- **200 rxs (25 ul)**

Highlights

The Master Mix mix can provide reliable results up to 50 real-time PCR cycles.

The master mix ensure the possibility of having up to three primer pairs in the same PCR reaction, with the respective 3 probes equipped with three different fluorophores.(for example:Taqman, MGB, LNA, Molecular Beacons, Scorpions)

Compatible instruments

| | |
|------------|--|
| No Rox | Bio-Rad iCyclers/ CFX96/ CFX 384, Roche Light Cyclers® QIAGEN/Corbett Systems Eppendorf Mastercylers®, Opticon 2 |
| Rox Dye I | Applied Biosystems 7000/7300/7700/7900, 7900HT, Applied Biosystems Step One Plus |
| Rox Dye II | Applied Biosystems 7500/ViiA7, QuantStudio 7 Flex ABI, Stratagene Real-time PCR Systems, Rotor-gene3000 |

Additional Material Required but not Supplied

1. Optical-grade qPCR tubes, plates, sealing films, and aerosol-resistant pipette tips
2. qPCR primers and probes
3. DNA or cDNA templates

Precautions

1. Fully thaw TaqMan Universal Multiplex qPCR Master Mix before use.
2. The TaqMan Universal Multiplex qPCR Master Mix contains glycerin. Mix gently before use without generating air bubbles. Spin briefly to collect all the contents at the bottom. After use, return it to -20°C immediately.
3. A Hot-start version of Taq polymerase is included in the master mix, allowing reaction setup at room temperature. After first thaw, the master mix is stable at 4 °C for 1 week.
4. Use the ROX reference dye according to the requirement of qPCR instrument to be used.
5. If applicable, use aerosol-resistant pipette tips to minimize contamination.
6. High quality DNA templates are recommended for optimal results.

Operation Description

Important points before reaction setup:

- (1) A final primer concentration of 0.2 μM is recommended for most reactions. However, to optimize individual reaction, a primer titration from 0.1 μM to 1.0 μM can be performed.
- (2) The length of amplified PCR products should ideally be in the range of 70-200 bp.
- (3) Prepare a serial dilution of the template to access standard curve and test primer efficiency.
- (4) Use 1 pg-50 ng of DNA template in a 20 μL reaction. The volume of template should not exceed 10% of the final PCR reaction volume.
- (5) Always include a no template control (NTC) reaction.
- (6) Triplicates are recommended as technical replicates in real-time PCR reactions.

Recommended Reaction

| Components | 25 μL |
|--|------------------------|
| TaqMan Probe 2X Multiplex qPCR Mater Mix | 12,5 μL |
| Forward Primer (10 μM) | 0.5 μL |
| Reverse Primer (10 μM) | 0.5 μL |
| Probe (10 μM) | 0.5 μL |
| 50X ROX Dye (as required by instrument guidelines) | 0.5 μL |
| DNA Template | 2,5 μL |
| Nuclease-free Water | Up to 25 μL |

Note :

- (1) Fully thaw the TaqMan Universal Multiplex qPCR Master Mix at room temperature, and gently mix well without creating bubbles. Spin down briefly in a microcentrifuge to collect all contents at the bottom.
- (2) Calculate the required volume of each components based on the number of reactions to be set up and add extra 10% volume of each component to compensate pipette errors.
- (3) Add all the common reaction components (primers and probes) in a master mix and mix thoroughly.
- (4) Dispense appropriate volumes of reaction mix into qPCR plates, and carefully seal it with an optical sealing film.
- (5) Add templates or NTC into wells containing the qPCR reaction mix.
- (6) Centrifuge the qPCR plates (tubes) at 2500 rpm to collect all the contents at the bottom of wells. The samples are ready for thermocycling.

Recommended PCR Program

| Step | Temp | Time | Cycles |
|-------------------------|------|-------|--------|
| Predenaturation | 95°C | 3 min | 1 |
| Denaturation | 95°C | 15 s | |
| Annealing and extension | 55°C | 30 s | 40-50 |

Note : (1) To ensure signal acquisition after extension, the extension temperature should be based on the T_m value of the primer probe. Line adjustment.

(2) It is recommended that the shortest predenaturation time should not be shorter than 3 min, and the longest should not exceed 10 min; the shortest denaturation time during the cyclic reaction is not less than 5s, and the longest is not more than 15s; the cyclic reaction; the shortest extension time in the application is not less than 10s, and the longest can be based on the primer probes and signals used by yourself. The set needs to be adjusted by itself.

Data Analysis

1. A standard curve is a linear regression analysis on the data plotted as the C_t values versus the log sample input concentration. If the standard curve correlation coefficient (R^2) is >0.98 , the template concentration data points are within the linear range of the assay. When the slope of a standard curve is between -3 and -3.5, the PCR amplification efficiency (E) is between 90 and 120%.
2. Ideally, the standard deviations (STD) of the C_t values between replicates should be <0.2 cycle. For most cases, the STDs of C_t values <0.5 cycle are acceptable.
3. The C_t value of a valid amplification should be less than the value of the NTC curve.

Storage: Upon receipt store at -20° for long term.

For Research Use Only

FISHER MOLECULAR BIOLOGY

36 Terry Drive Trevose,
PA 19048 - USA