

PLASMID DNA ENDOFREE EXTRACTION MIDI PREP KIT DE-50103


The Plasmid DNA Endofree Extraction Midi Kit is designed for rapid and efficient extraction of high quality plasmid DNA from Bacterial cultures with Spin Columns. This kit is suitable for extracting plasmid DNA from 5-12 mL bacterial cultures. After the bacterial cells are treated with alkaline lysis, high salt, and low pH conditions, the plasmid is released and specifically and efficiently adsorbed by the silica gel membrane of the centrifugal column. The bacterial endotoxins can be removed by washing with Buffer ToxOut. Then, proteins and other impurities can be removed by washing with Buffer W2. Finally, high - purity endotoxin - free plasmid DNA is eluted under low - salt and high - pH conditions. The purified plasmid DNA is suitable for various routine operations, including conventional cell transfection, restriction enzyme digestion, PCR, sequencing, ligation, transformation, library screening, in vitro Translation, etc.

Specifications:

High yield: Up to 70 µg plasmid DNA can be quickly extracted from 5 mL-12mL bacterial cultures within 1 h.

High purity: Optimized Buffer is used to reduce endotoxin, protein, SDS and other residues.

Wide application: Plasmid DNA prepared is suitable for a variety of downstream applications including restriction enzyme digestion, PCR, sequencing, ligation, transformation and cell transfection.

Components	DE-50103 50 preps	 FMB vs Competitor Plasmid DNA extracted from Bacterial Culture cells
BL2 Buffer	25 ml	
S1 Buffer	30 ml	
S2 Buffer	30 ml	
S3 Buffer	30 ml	
ToxOut Buffer	30 ml	
W2 Buffer	15 ml	
EB2 Buffer	15 ml	
RNase A* (10 mg/mL)	300ml	
Spin Column 5	50 pk	
Collection tubes (2ml)	50 pk	

Precautions:

1. RNase A is provided in the kit. After brief centrifugation, add all the RNase A to Buffer S1, mix well, and store at 2 - 8 °C. Before use, let it stand at room temperature for a while until it returns to room temperature.
 2. Before the first use, add 45 mL of absolute ethanol to Buffer W2. After use, immediately tighten the bottle cap to prevent the evaporation of ethanol.
 3. Each time before use, check whether precipitates have formed in Buffer S2 and Buffer S3. If there are precipitates, dissolve them at 37°C before use.
 4. The general incubation time for bacteria is 12 - 16 hours. If the inoculum size is large, the incubation time should be reduced.
- Over - cultivation may reduce the quality of the plasmid and even lead to plasmid DNA mutation.
5. The yield of the plasmid is closely related to the amount of bacteria, plasmid copy number, plasmid size, and the degree of standardized operation.
 6. If the plasmid to be extracted is a low - copy plasmid, a large plasmid larger than 10 kb, an Agrobacterium plasmid, or a Gram - positive bacterium plasmid, increase the amount of bacterial cells. Use 5 - 12 mL of an overnight culture. At the same time, increase the amounts of Buffer S1, Buffer S2, and Buffer S3 proportionally. Pre - heat the elution buffer in a 60°C water bath. Appropriately extend the adsorption and elution time to increase the extraction efficiency. The other steps remain the same.

Cat.#	Description	Size
DE-50103	Plasmid DNA Endofree Extraction Midi Prep	50 preps

Storage: The kit is stable for 12 months when store at RT

For Research Use Only

FISHER MOLECULAR BIOLOGY

36 Terry Drive
Trevose, PA 19048 - USA

General Protocol: Operational instructions

1. Column equilibration: Add 400 μ L of Buffer BL2 to Spin Columns 5 (place the adsorption column in the collection tube), centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 1 min, discard the filtrate in the collection tube, and put the adsorption column back into the collection tube.

Note: Try to use the column treated with the eBuffer BL2 on the same day to avoid affecting the subsequent experimental results due to long - term storage.

2. Take 5 - 12 mL of an overnight - cultured bacterial liquid, centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 1 min at room temperature, and remove the supernatant as completely as possible.

Note: If there is a large amount of bacterial liquid, it can be added to the tube in batches for centrifugation to collect the bacterial cells. Do not use too much bacterial cells, otherwise the lysis will be insufficient.

3. Add 500 μ L of Buffer S1 (please check if RNase A has been added first), and then vortex or pipette thoroughly to resuspend the bacterial pellet evenly.

Note: The bacterial cell pellet must be suspended evenly. If there are incompletely suspended bacterial lumps, it will affect the lysis, resulting in a decrease in the concentration and purity of the extracted plasmid.

4. Add 500 μ L of Buffer S2, gently invert and mix to completely lyse the bacterial cells until the solution becomes clear and viscous.

Note: Invert and mix gently, do not shake vigorously to avoid breaking the genomic DNA. At this time, the bacterial liquid should become clear and viscous, and the time used should not exceed 5 min to prevent the plasmid from being damaged. If the solution does not become clear, it may indicate that the amount of bacteria is too large. You can increase the amount of Buffer S2, and the amount of Buffer S3 should also be increased accordingly in the subsequent operations.

5. Add 500 μ L of **Buffer S3**, and immediately invert gently to mix. A red-yellow precipitate will form. Continue mixing until the solution turns completely yellow. Let it stand at room temperature for 2 minutes, then centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 5 minutes.

Note: After adding Buffer S3, mix immediately to avoid local precipitation. If there are still purple floating substances in the supernatant, it means that the renaturation is insufficient. Continue mixing until the solution completely turns into a clear yellow color.

6. Carefully transfer the supernatant to a new centrifuge tube. Add 0.3 times the volume of isopropanol and mix. Then transfer the solution to the adsorption column and let it stand for 2 minutes to allow the plasmid DNA to bind fully with the silica membrane in the column. Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 s, and discard the filtrate in the collection tube.

Note: The maximum volume of the adsorption column is 750 μ L. The remaining liquid should be passed through the column in batches.

7. Add 500 μ L of **ToxOut Buffer** to the adsorption column, let it stand at room temperature for 5 min, and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 1 min. Discard the filtrate in the collection tube.

8. Add 600 μ L of Buffer W2 (please check if absolute ethanol has been added first), centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 second at room temperature, and discard the filtrate in the collection tube.

Note: Buffer W2 is a concentrated solution. Add absolute ethanol as required. After use, immediately tighten the bottle cap to prevent the evaporation of alcohol.

9. Add 500 μ L of Buffer W2, and centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min at room temperature to remove the residual liquid.

Note: This step cannot be omitted, otherwise the residual ethanol will affect the downstream experiments.

10. Place the adsorption column in a new 1.5 mL centrifuge tube (prepared by yourself), add 100 - 300 μ L of elution buffer

Buffer EB2 to the center of the adsorption membrane, let it stand at room temperature for 2 min, and centrifuge at 12,000

rpm ($\sim 13,400 \times g$) for 1 min to collect the plasmid solution in the centrifuge tube.

Note: To increase the elution efficiency, Buffer EB2 can be pre - heated at 60 ° C. Buffer EB has a single component and will not affect downstream restriction enzyme digestion, transfection, or other molecular biology experiments. Please use it with confidence. Elution with deionized water for elution, adjust its pH value to between 8.0-8.5 with NaOH before use. To increase the plasmid recovery rate the obtained solution can be added back to the centrifuge tube, let it stand at room temperature for 2 minutes and centrifuge again for collection.

11. The obtained DNA should be stored at -20°C or used directly for subsequent experiments.

For Research Use Only

FISHER MOLECULAR BIOLOGY

36 Terry Drive
Trevose, PA 19048 - USA