

Taqman Probe qPCR

TAQMAN PROBE FAST Qpcr MASTER MIX


The TaqMan Probe Fast qPCR Master Mix is an optimized, ready-to-use PCR mixture of Hot-start Taq DNA Polymerase, PCR buffer, Magnesium and dNTPs, except DNA template and primers.

Compatibility

No ROX	ABI 7000, 7300, 7700, 7900HT and 7900HT Fast, StepOne, StepOne Plus
Low ROX	ABI 7500, 7500 Fast, ViiA 7, QuantStudio; Roche LightCycler; Stratagene Mx3000, Mx3005P and Mx4000
High ROX	BioRad: iCycler, MyiQ, MiQ 2, iQ5, CFX-96, CFX-384, MJ Opticon, Option2, Chromo4, MiniOpticon Qiagen: Roto-Gene Q, Roto-Gene3000, Roto-Gene 6000 Eppendorf: Mastercycle realplex - Illumina: Eco RealTime PCR System Cepheid: SmartCyler- Roche:LightCycler 480,LightCycler 2.0

Cat.No.	Size
FS-T-1072F	5 X 1 ml
ROX Dye (1x)	1 vial

1ml = 100 Reactions



DATASHEET

Per ordini e/o offerte contattare: Società Italiana Chimici

Tel: +39 06 8818936 – 8800211 Email: offerte@sichim.com

Taqman Probe qPCR

TAQMAN PROBE qPCR UDG

Entrans 2X qPCR Probe Master Mix with UDG is a ready-to-use reagent for probe-based qPCR reactions, containing all components except primers, probes and templates.

Compatibility

No ROX	ABI 7000, 7300, 7700, 7900HT and 7900HT Fast, StepOne, StepOne Plus
Low ROX	ABI 7500, 7500 Fast, ViiA 7, QuantStudio; Qiagen: Roto-Gene Q, Roto-Gene3000, Roto-Gene 6000
High ROX	BioRad: iCycler, MyiQ, MiQ 2, iQ5, CFX-96, CFX-384, MJ Opticon, option2, Chromo4, MiniOpticon Roche: LightCycler 480, LightCycler 2.0 Eppendorf: Mastercycle realplex - Illumina: Eco RealTime PCR System Cepheid: SmartCycler

Contents	Size	Price
FS-T-50217	500 RXN	
Rox Dye I (high Rox) 50X	1 vial	CC
Rox Dye II (high Rox) 50X	1 vial	

1ml = 100 Reactions

DATASHEET →

